



En

User Manual

VISO10 – MSF (United Kingdom)


Please refer to our website: **www.geemarc.com** for an up-to-date user guide - there may be important updates and changes you need to be aware of.

For best use of the Viso10 – please read all sections of this user guide. Keep the user guide in a safe place for future reference.

1. Features

- Easy to read and understand time display
- Calendar up to year 2099
- Day of the week & month in 8 languages (user selectable)

2. Unit Appearance



Display (LCD)	Push-Button Controls
A1: Day of the Week A2: Radio Control Icon A3: Time A4: Month A5: Date A6: Year	B1:"MODE/ Button B2: "▼" Button B3: "▲" Button B4: "WAVE/ Button (Restart the signal search)
A7: Locked keys A8: AM / PM A9: Season indicator Winter / Summer	
Structure	
C1: Hanging Keyhole	

3. Getting started:

- Open the battery compartment cover (C2).
- Insert 3 x AA alkaline batteries. Observe the correct polarity "+" and "-". A "click" is heard when the batteries are inserted correctly.
- Replace the battery compartment cover.

4. Time Setting

4.1 Manual Time Setting:

- In time display, press and hold "MODE/ button(B1) for 3 seconds to enter Clock/Calendar setting Mode.
- Press "▼" (B2) or "▲" (B3) button to adjust the setting then press "MODE" button (B1) to confirm each setting.
- (Press and hold the "▼" (B2) or "▲" (B3) buttons to reach the desired value faster)
- The setting sequence is shown as follows: Weekday & Month Language, 12/24 Hr, RCC (Radio Clock) ON/OFF, Time Zone, Hour, Minute, Year, Month, Day, DST.
- There are 8 languages for day of the week and month: English, German, French, Spanish, Italian, Dutch, Polish and Swedish.

4.2 Radio Controlled Time:

• When Viso10 is switched on, it will automatically scan the MSF time signal. The Radio Control Icon "((())" (A2) flashes on the (LCD) display – as indicated below.

(A2) flashes	Receiving MSF signal
(A2) ON	MSF signal received successfully
(A2)disappears	MSF signal reception failed

• To maintain accurate time, the clock will

C2: Battery Compartment Cover

C3: Stand

automatically synchronize with the MSF radio signal every day at 2.00 am. If synchronization with the MSF radio signal fails, the "(())" symbol on the display will disappear. The clock will then attempt to synchronize with the MSF radio signal at 3am and 4am.

- To scan the MSF radio signal manually, press and hold the "((())" button (B4) for 3 seconds. Viso10 scans for reception of the MSF signal every 7 minutes. If reception fails, the scanning will stop. ("((()))" on LCD disappear). Scanning will repeat again on next full hour. Example: if scanning failed at 8:20 am, it will scan again at 9 am.
- To stop scanning MSF radio signal, press and hold the "(")" button (B4) for 3 seconds.
- "DST" will show on the display to indicate time mode in winter time.

NOTE:

(1) The clock will automatically exit from all setting modes after approx. 15 seconds without any adjustment.

5. Lock buttons

Press and hold "▼" (B2) and "▲" (B3) buttons at the same time for 5 seconds to lock all other button controls. The 'padlock' icon "₁" is displayed to show that the button controls are locked. Press and hold "▼" (B2) and "▲" (B3) buttons again for 5 seconds to unlock.


6. Precautions

- The clock will lose its time information when the battery is removed.
- Do not expose the unit to direct sunlight, extreme heat or cold, high humidity or wet

areas.

 Never clean the device using abrasive or corrosive materials or products. Abrasive cleaning agents may scratch plastic parts and corrode electronic circuits.

7. Radio controlled time reception range: up to 1200km from MSF transmitter.



8. Guarantee

From the moment you purchase the Geemarc Viso10V2, Geemarc will guarantee it for the period of two years. During this time, all repairs or replacements (at our option) are free of charge. Should you experience a problem then contact our helpline or visit our website at www.geemarc.com. The guarantee does not cover accidents, negligence or breakages to any parts. The product must not be tampered with or taken apart by anyone who is not an authorised Geemarc representative. The Geemarc guarantee in no way limits your legal rights.

IMPORTANT: YOUR RECEIPT IS PART OF YOUR GUARANTEE AND MUST BE RETAINED AND PRODUCED IN THE EVENT OF A WARRANTY CLAIM.

Declaration: The declaration of conformity may be consulted at **www.geemarc.com**

RF power	< 250mW
Power requirements	AA batteries (3)

9. Recycling directives

The WEEE (Waste Electrical and Electronic Equipment) has been put in place for the products in end of life are recycled in the best way.

When this product is end of life / no longer working, please do not put it in your domestic waste bin.

Please use one of the following disposal options:

- Remove the batteries and deposit them in an appropriate WEEE skip.
- Deposit the product in in a appropriate WEEE skip.
- Or hand the old product to the retailer. If you purchase a new one, they should accept it.



For product support and help visit our website at www.geemarc.com

E-mail : help@geemarc.com telephone 01707 387602

lines are open open 09.00 to 16.00 Mon to Fri Parc de l'étoile, 2 rue Galilée 59760 Grande-Synthe, France

Made for Geemarc Telecom S.A. in China



UGViso10V2_En_A4_V1