


How to use your Kylie® bed pad:

- Place your Kylie® bed pad directly on the made up bed on top of the bedsheet.
- Kylie® Pad should be placed with the Blue/Pink soft quilted fabric facing up.
- Bottom non slip waterproof barrier will grip to your sheet holding your pad firmly in place.
- Tuck the wings of your Kylie® bed pad firmly under the mattress.









Composition

- Soft quilted 100% polyester acquisition layer.
- Highly absorbent soaker layers, polyester / viscose.
- Waterproof polyurethane membrane.
- Polyester backsheet.
- Fitted with 50cm polycotton wings / flaps.

Kylie® bed pad washing instructions

Kylie® Bed Pad must be washed prior to use, this will ensure that absorbency and comfort are at their optimum.

-  Do not use fabric conditioner
-  Do not iron
-  Do not dry clean
-  Machine wash
-  Do not bleach
-  May be tumble dried low/medium



Capatex Code	Colour	Absorbency	Barcode	Product Weight	Pack Dimensions
(with wings) IN-BP-KYLIE1-PINK-W	■	1 Litre	5060077030190	0.35Kg	240 x 300 x 50 cm
(with wings) IN-BP-KYLIE2-PINK-W	■	2 Litres	5060077030206	0.65 Kg	240 x 300 x 70 cm
(with wings) IN-BP-KYLIE3-PINK-W	■	3 Litres	5060077030213	0.80 Kg	240 x 300 x 85 cm
(with wings) IN-BP-KYLIE4-PINK-W	■	4 Litres	5060077030237	1.15 Kg	220 x 310 x 110 cm
(with wings) IN-BP-KYLIE5-PINK-W	■	5 Litres	5060077032880	1.25 Kg	220 x 310 x 120 cm
(with wings) IN-BP-KYLIE1-BLUE-W	■	1 Litres	5060077031425	0.35Kg	240 x 300 x 50 cm
(with wings) IN-BP-KYLIE2-BLUE-W	■	2 Litres	5060077031432	0.65 Kg	240 x 300 x 70 cm
(with wings) IN-BP-KYLIE3-BLUE-W	■	3 Litres	5060077031494	0.80 Kg	240 x 300 x 85 cm
(with wings) IN-BP-KYLIE4-BLUE-W	■	4 Litres	5060077031500	1.15 Kg	220 x 310 x 110 cm
(with wings) IN-BP-KYLIE5-BLUE-W	■	5 Litres	5060077032897	1.25 Kg	220 x 310 x 120 cm