



TL 100 Tageslichtlampe 2in1

Licht: Wohlbefinden für Körper und Seele

Zur Anwendung bei Lichtmangelerscheinung in den Wintermonaten: Unausgeglichenheit, gedrückte Stimmung, Energie- und Antriebslosigkeit

2 in 1: Tageslichtlampe und Stimmungslicht

Komfortabel und einfach über „beurer LightUp“ App bedienbar

Stimmungslicht mit Farbwechsel und 256 Farben

Schlankes Design dank LED -Technologie

Dimmer mit 10 Helligkeitsstufen und Memory-Funktion

Simulation von Sonnenlicht:
Lichtstärke ca. 10.000 Lux (12 cm)

LED-Timer zur Anzeige der Behandlungsdauer

Behandlungszeit-Anzeige in 15-Min-Schritten

Zertifiziertes Medizinprodukt (Tageslichtlampe)

Besonders helle und gleichmäßige Ausleuchtung

Beleuchtungsfläche ca. 33 cm Ø; 855 cm²

Energiesparende LED Technologie

Flackerfreier Start dank LED-Technologie -
sofort angenehmes Licht ohne Flimmereffekte

Spektrum: 400-800 nm

Farbwiedergabeindex: > 90

Farbtemperatur: 5300 K +/- 300 K

Material:

Leuchtschirm: PC

Gehäuse und Standfuß: ABS

Eingang: 100-240V~, 50-60Hz, 1,5A

Ausgang: 24V --- 1,5A

Produktmaße (LxBxH): 306 x 335 x 203 mm

Produktgewicht: ca. 1300 g, Netzteil ca. 173 g

3 Jahre Garantie

VE: 2 / Transportkarton: -

EAN-Nr.: 4211125 60835 4

Art.-Nr.: 608.35

*Basisfunktionen direkt am Gerät einstellbar.

TL 100 Brightlight 2in1

Light: well-being for body and soul

For use in case of light deprivation symptoms in the winter months: e.g. low mood, lack of energy and drive

2-in-1: brightlight and mood light

Convenient and easy to use via "beurer LightUp" app

Mood light with colour change and 256 colours

Slim design thanks to LED technology

Dimmer with 10 brightness levels and memory function

Simulation of sunlight:
Light intensity of approx. 10,000 lux (12 cm)

LED timer to display treatment time

Treatment time display in 15 minute intervals

Certified medical product (brightlight)

Exceptionally bright and even illumination

Illumination surface approx. 33 cm Ø; 855 cm²

Energy-saving LED technology

Flicker-free start-up thanks to LED technology -
immediately produces pleasant light without flickering

Spektrum: 400-800 nm

Color Rendering Index (CRI): > 90

Color temperature: 5300 K +/- 300 K

Material:

Lamp shade: PC

Housing and stand: ABS

Input: 100-240V~, 50-60Hz, 1.5A

Output: 24V --- 1.5A

Product measurements (LxWxH): 306 x 335 x 203 mm

Product weight: approx. 1300 g, mains part approx. 173 g

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 60835 4

Item no.: 608.35