



110

The 110's slim design makes it the perfect solution for older or narrow staircase evacuations. The extension handle can be adjusted to ensure a comfortable operating position and the footrest provides a safe place for the passenger's feet. The unique friction belt system provides a smooth evacuation for the mobility impaired.

FEATURES:

- 1 Multi height extension handle**
For ease of horizontal movement (use fully extended for descent).

- 2 Adjustable passenger head strap**
Limits movement, securing the passenger in place.

- 3 Sewn-in safety belt**
Improves safety and durability.

- 4 Unique friction belt system**
Self-braking track to control smoothness and speed of descent.

- 5 Redesigned kickstand, with repositioned rear wheels**
Increased stability and enhanced operator handling.

- 6 Large 100mm castor wheels**
Ideal for operation over uneven terrain.

- 7 Foldable Footrest**
For added passenger comfort and safety.

SPECIFICATION:

Height	1020mm (folded)	1305mm (unfolded)
Width	460mm (folded)	460mm (unfolded)
Depth	245mm (folded)	1075mm (unfolded)
Weight (unladen)	11.5kg	-
Payload capacity	150kg	-
Stair angle	28° - 40°	-



APPROVALS/ACCREDITATIONS:

- + TUV NORD
- + CE & UKCA Mark
- + Registered Class 1 Medical Device
- + ISO 13485: Medical Device Quality Management
- + ISO 9001: Quality Management System

BENEFITS:

- + Easy to use
- + Narrow staircase operation
- + One-person operation
- + 23 point quality check throughout manufacture
- + Lightweight

INCLUDES:

- + Photoluminescent signage
- + Wall mounted hooks
- + Dust cover
- + User guide

