

# EV4000

## The standard model evacuation chair

Provides a quick, easy and safe emergency evacuation with a 200kg (440lbs) weight capacity and allows for single person operation.



### FEATURES

- 1. Adjustable self braking tracks**  
Reinforced adjustable tracks control the speed of descent
- 2. Anti-slip height adjustable handle**  
For ease of deployment, operation and compact storage
- 3. Head support with immobilisation strap**  
Limits movement, securing the passenger
- 4. Quick release 4-point harness**  
4 point harness for increased passenger welfare
- 5. Heavy duty front wheels**  
For smoother operation over rugged terrain
- 6. Industrial locking castor wheels**  
For situations where the chair needs to be stationary
- 7. Way finding decals**  
For easy location in emergency situations
- 8. Adjustable seating position**  
Allows for the ideal passenger seating position, either flat seat or hammock

### DIMENSIONS

Height	Width	Depth	Weight	Weight Capacity
124cms 48"	57cms 22.5"	26.5cms 10.5"	14.3kgs 31.5lbs	200kgs 441lbs

### BENEFITS

- Single person operation up to 200kg (441lbs)
- Easy to open, ready to operate in seconds
- Provides a rapid solution for safe emergency evacuation
- Simple to operate
- Lightweight construction
- Minimal moving and handling required
- Allows small operators to transport passengers larger than their size
- Foldable for easy storage
- Fire retardant material
- Wall mountable
- Low maintenance

### INCLUDES

- LIFETIME WARRANTY
- 4 Point Passenger Harness
- Cushioned Head Support
- Adjustable Padded Seat
- Way Finding Signage
- Wall Mounting Hooks
- Protective Cover
- User Guide • FREE UK Delivery