

# EV9000

3 speed powered stair climber with increased weight capacity of 227kgs for both upward and downward evacuation

Our powered solution uses Power-Climb technology with an increased weight makes negotiating up and down stairwells that bit more less stressful. Single person operation.



## FEATURES

- 1. Anti-slip height adjustable handle with head cushion**  
For ease of operation and compact storage, cushion adds comfort for passenger during transit
- 2. Easy to use 3 speed control panel with safety light**  
With on/off power and ascent and descent buttons
- 3. Additional operator handles**  
For added stability when employing additional operators
- 4. Quick release 4-point harness**  
4 point harness for increased passenger welfare
- 5. Rechargeable lithium ion battery pack**  
Charging time 8 hours (When fully charged is the equivalent to going up and down a 5 storey building flight of stairs 30 times)
- 6. High strength aluminium frame**  
Increased strength and robustness for continued use
- 7. Adjustable self braking tracks**  
Reinforced adjustable tracks control speed of descent on all types of stairways
- 8. Industrial locking rear wheels**  
For application in situations where the chair needs stay in a stationary position

## DIMENSIONS

Height	Width	Depth	Weight	Weight Capacity
116cms	57cms	32cms	33kgs	227kgs
45"	22.4"	12.6"	72.7lbs	500lbs

## BENEFITS

- Increased weight carrying capacity up to 227kg (500lbs)
- Easy to open, ready to operate in seconds
- Simple to operate
- Single person operation
- Lightweight construction
- Minimal moving and handling required
- Allows small operators to transport passengers larger than their size
- Foldable for easy storage
- Fire retardant material
- Wall mountable
- Non sparking electrical components eliminating risk of fire
- High capacity battery with burnout protection

## INCLUDES

- Battery Charger (UK/EU/US)
- Rechargeable Battery
- 2 Keys • User Guide
- Passenger Headrest
- 4 Point Passenger Harness
- Way Finding Wall Sign
- Wall Mounting Hooks • Protective Cover
- Manufacturers Warranty