

## COMPACT EVACUATION MAT



### PRODUCT CODE: EVM-W

Our new Evacuation Mat is designed to be small and compact as well as being very hard wearing. The unit is constructed of highly durable fabric with a foam pad for comfort. The principles of the unit are the same as all our evacuation products. The webbing runs the full length of the unit encapsulating both head and foot end pulling straps in one piece of very strong webbing. Together with the two body straps for tall and short patients along with a leg strap. Once the patient is on the unit they can be enveloped for security by the wings on either side. Secured by the webbing and buckles and then pulled to safety by means of the head and foot pull straps.