

sharpsafe®

The all new 5L Sharps Container

With DUAL PROTECT **PRO™**

DATA SHEET



The Sharpsafe® 5L container is ideally suited for sharps disposal in medium use areas and ideal for acute settings. Offering industry leading features such as the new Dual Protect Pro™, the container is designed to encourage horizontal fill to maximise capacity. The new container design has a 13% greater capacity than the previous model. Combined with the low carbon footprint recycled material in the base, this container supports carbon reduction goals within healthcare settings.

Lid colour options:

Yellow

Product Code: 51025430 GTIN Code: 5056271508955

Orange

Product Code: 51025410 GTIN Code: 5056271508962

Purple

Product Code: 51025420 GTIN Code: 5056271508979

Blue

Product Code: 51025440 GTIN Code: 5056271509327

Blue Scottish

Product Code: 51025220 GTIN Code: 5056271509884

Red Scottish

Product Code: 51025280 GTIN Code: 5056271509891

Key Features:

- ✓ The first sharps container to feature Dual Protect Pro™ 1
- ✓ In-Mould Label
- ✓ 100% Recycled Plastic Base
- ✓ Luer Lock/Luer Slip Disconnection Feature 2
- ✓ Star Disconnection Feature 3
- ✓ Integrated Handle
- ✓ Customisable Label (on request)
- ✓ Four Click Assembly

Safer today for a sustainable tomorrow.

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. © Copyright Vernacare International Limited 2024. All rights reserved. CWM008 REV02



Vernacare
clinical waste
management

Pocket
Small
Medium
Large
Specialist
Accessories

The all new 5L Sharps Container

Certifications	
ISO 23907-1:2019	Yes
BSI Kitemark	Yes
UN3291	Yes
ISO 9001	Yes
ISO 14001	Yes
ISO 45001	Yes
Medical Device	No
CE Marked	N/A

Conformities	
AS 23907	Pending
TRBA 250	Yes

Logistical Data	
Carton Qty	48
Weight When Empty (kg)	0.82
Net Weight (kg)	16.7
Carton Width (mm)	399
Carton Height (mm)	490
Carton Length (mm)	593
Volume (m ³)	0.1159
MOQ	1 Carton
Lead Time	10 Days
Storage Temperature Conditions	Ambient
UK Pallet – Cartons Per Pallet	20
UK Pallet – Units Per Pallet	960
Euro Pallet – Cartons Per Pallet	16
Euro Pallet – Units Per Pallet	768

Product Information	
Total Volume (L)	5.19
Usable Volume (L)	4.31
Capacity (%)	83
Weight of Base (g)	217.7
Weight of Lid (g)	113
Combined Weight (g)	330.7
Fill Line Height (mm)	211.6
Height (mm)	260.6
Width (mm)	170
Depth (mm)	205
Aperture Size (mm) (L x W)	61 x 65
Sterile	No
Lot/Batched	Yes
Expiry Date	N/A

Material Information	
Base Material	Recycled Polypropylene
Lid Material	Polypropylene
Latex Free	Yes
Free From Animal Products	Yes
REACH Statement	Yes

Accessory Compatibility	
Accessory	Product Code
Stand	42570000
SafeClip Bracket	45000000
Wall Loop	Coming Soon

Pocket
Small
Medium
Large
Specialist
Accessories

Safer today for a sustainable tomorrow.

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. © Copyright Vernacare International Limited 2024. All rights reserved. CWM008 REV02



Vernacare
clinical waste
management